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Advan pursues an innovative vision of implantology: the implant system is placed at the center of the digital workflow as an essential requirement for any dental team in order to meet at best the clinical needs of patients.

With this purpose, **Advan** is the result of the union of two companies with more than 25 years of experience in design, production and sale of dental implant products and with the aim of distributing their products in Italy and worldwide.

Advan can fully follow the development of a medical product, from concept to engineering, from mechanical tests to management of clinical trials up to certification and mass production.

This certification is issued for the total quality system, all activities are therefore according to the operating procedures of the quality manual, the production processes constantly monitored to ensure the quality and safety of the products and their traceability.

Operational activities are carried out in two separate production sites: precision machining at the plant, exclusively dedicated to implantable medical products, equipped with the most up-to-date CNC machines for providing a product of superior precision, strictly duplicated for every lot and controlled with the most advanced systems of metrology (laser technology); surface treatments, assembly of semi-finished products, marking, decontamination and packaging are carried out at the production site equipped with a certified cleanroom for the production of implantable medical products (ISO 7 Class Clean Room).

Advan also provides its know-how and operative structures for companies or professionals who want to carry out their projects, offering all the necessary resources and support from the development of prototypes to production.

In **Advan** the experience is combined with the constant and meticulous attention concerning the materials development, the passion for innovation, the most up-to-date production technologies and the clinical guidelines, supported by the scientific contribution of our Opinion Leaders.

Mission

To study and to generate implants systems that are easy to use, with high clinical performances and highly predictable long term results, also through the digital workflow, are the clear and defined objectives that everyone, in **Advan**, knows and constantly pursues.

Advan is committed to deliver significant value to Patients and to Professionals, everyone is personally responsible and aware that quality of products touches the quality of life for thousands of lives every day.

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The **ONE Implant System** has been designed with the over 25 years experience of **Advan** in the design, production and sale of dental implants and instruments dedicated to implant prosthetic rehabilitation.

ONE is synonymous with accuracy, simplicity, compatibility and versatility contained in a single dental implant.

The **ONE implant** screw has a modern and innovative design and, thanks to its peculiar apical treads, guarantees an easy positioning and an excellent primary stability even if it's a very user friendly implant system, but performing with advanced surgical techniques too (i.e. reaching bicorticalism in the immediate loading of the upper jaw or finding great primary stability in immediate post-extraction sockets). The atraumatic apex of the **ONE** implant screws allows mini-elevation of the maxillary sinus even without the aid of bone grafts. **Advan's** experience with high-end implant systems has allowed to use innovative features even in a more easy to use implant system: the crestal module (to respect the marginal bone) and the fluid discharge channel (to respect the alveolar bone) are two significant examples of the know-how transfer, making these specifications available to all professionals.

The **ONE system** is characterized by a single surgical kit and a single fixture design but available with three different connections:

- **Internal hexagonal**, to easily manage most of the clinical cases, from single implant to full arch rehabilitation, without changing the professional's habits (connection compatible with the most widespread internal connection system in the world);
- **Internal conical**, to better manage the hard and soft tissues in the aesthetic sectors or when the tissue availability requires an accurate management of the connective gingival (connection compatible with the **Advan GTB Implant System**);
- **External hexagonal**, for money saving implant rehabilitations such as Toronto Bridge or Overdenture (connection compatible with the most widespread external connection system in the world), available in 2020.

The drills sequence is easy and intuitive, does not require specific rotating instruments and is based on the most widespread drills sequence so that **ONE Implant System** can be easily used even by professionals who has recently invested in a new surgical kit or that has just replaced the drills.

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
The prosthetic connection, regardless of the type of connection chosen, guarantees a native platform shifting present on all implants diameters. Furthermore, the single connection diameter makes the implant prosthodontics quick and save.

The prosthetic components have been carefully selected to allow you to work with ease but maintaining a high aesthetic&funcional standard, that's all is needed by modern dental practice. The dental practice should always choose an implant systems that allow to achieve the pre-established result using just a limited number of prosthetic components, everything perfectly integrated into the most widespread digital diagnostic systems and digital prosthodontics production: **ONE is a Native Digital System.**

Advan's experience in managing the digital workflow allows all our customers to update their workflow or to ask for a perfect integration of Advan's Implant System with the workflow already established.

The great versatility and compatibility makes **ONE Implant System** ideal for professionals who need to switch the implant prosthetic connection in relation of the prosthetic rehabilitation. ONE Implant System becomes an incredibly efficient and effective tool for dental practice with multiple implant surgeons, each with its own philosophy of work, ensuring a standardization of workflow and material supply.

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
BIBLIOGRAFIA

BIBLIOGRAPHY

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- (1) Baggi L, Di Girolamo M, Mirisola C, Calcaterra R: Implant Dent. 2013 Aug;22(4):344-50. doi: 10.1097/ID.0b013e3182943062. Microbiological evaluation of bacterial and mycotic seal in implant systems with different implant-abutment interfaces and closing torque values.
- (2) Baggi L, Cappelloni I, Di Girolamo M, Maceri F, Vairo G: J Prosthet Dent. 2008 Dec;100(6):422-31. doi: 10.1016/S0022-3913(08)60259-0. The influence of implant diameter and length on stress distribution of osseointegrated implants related to crestal bone geometry: a three-dimensional finite element analysis.
- (3) Rompen et al.: J Prosthet Dent. 2007 Jun; 97(6 Suppl): S119-25
Soft tissue stability at the facial aspect of gingivally converging abutments in the esthetic zone: a pilot clinical study
- (4) Guo X E 2001
Biomechanics and Remodeling of Trabecular Bone, in Recent Advances in Biomechanics, Edited by C. Zhu and M. Long, Chinese Higher Education Press and Springer-Verlag, pp 120-129
- (5) Biasotto M, Di Lenarda R, Cadenaro M:
Effects on titanium surface topography of different surface treatments; J Dent. Res. 2000, 79 (special issue), abstr. 1757
- (6) Sandrucci M A, Basa M, Biasotto M, Bareggi R:
Fibronectin (FN) arrangement and collagen expression in human fibroblast cultured on titanium surface; Italian Journal of anatomy and embryology (Liv congress società di anatomia, Firenze 21-14 settembre 2000), 2000 vol.105, suppl.1:209

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THE ONE DENTAL IMPLANT
IS AN OSTEOINTEGRATED
SYSTEM WITH CONTROLLED
AND GRADUAL
TRANSMISSION OF THE
MASTICATORY LOAD
SUITABLE FOR IMMEDIATE
LOAD PROTOCOLS

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Cylindrical profile of the implant neck
Minimizes the stress of the cortical bone

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Fluids discharge channel

*Reduction of the hydrostatic compression
of the bone clot*

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Conical profile of the implant root

*Ideal distribution of the masticatory
load, minimal tensile and maximum
compressive forces on the bone*

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Double, progressive and modified reverse buttress main thread

*Faster insertion, better primary stability,
ideal load transmission. Enhanced BIC*

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
Double apical thread perfectly filled with main tread

*Faster insertion, better primary stability,
under preparation in soft bone, no
compression in compact bone*

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The implant screws are processed, cleaned and treated accordingly to the **certified procedures**. The final decontamination is performed by **cold plasma of argon** inside the **Clean Room** and then the implant screws are immediately stored.

The implant screws are removed from storage only inside the **Clean Room** of the production factory. In this way, the superficial oxide layer will be formed in a **controlled atmosphere**, not allowing the slightest contamination of the implant surface.

After packaging in the **Clean Room**, the implant screws are sent to the **sterilization process by β rays**.

The packaging is characterized by:

- **Decontamination** by cold plasma of argon
- **Sterilization** by β rays.
- **Double sterility**: outer barrier provided by blister PET-G with sheet in TYVEK®; internal barrier provided by the vial in polycarbonate
- **Cover screw** positioned in the vial cap under a seal in TYVEK®
- **Shaped blister** to avoid impacts to the vial
- **Implant screws in contact only with elements in titanium**, to avoid phenomena of bimetalism or contamination from contact with plastic material

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


All the informations on the implant blister package label are replicated on the adhesive labels for the patient record.

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
Choose the implant type, length and diameter and take the blister out of the cardboard box.

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
The vial containing the implant is sterile and lodged in the blister. The product description and the lot number are indicated on the label. Open the blister.

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

Surgical cover screw is placed in the vial cap and sealed with a Tyvek film.

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
Gently open the vial cap (do not pull up with force).

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
Connect the Implant Driver (choose the implant driver matching with the implant connection) and screw at low speed (10-15 rpm).

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
Pick the cover screw from the vial cap and screw it onto the implant with the Prosthetic Driver (07-EG05; 07-EG10; 07-EG20; 02EGM10).

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
SURGICAL KIT ONE

The surgical kit consists of a single box containing all the surgical and prosthetic instruments necessary for the realization of the osteotomy, the implant positioning and to its prosthesis. Including Radel made surgical tray (non-toxic, autoclavable, metal-free). The surgical kit have to be paired with the Implant Driver Set corresponding to the implant connection chosen.

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
ONE INTERNAL IMPLANT DRIVER SET

07NDD02
07NDM10
07NDM11

ONE INTERNAL DRIVER HANDPIECE LONG
ONE INTERNAL DRIVER MANUAL SHORT
ONE INTERNAL DRIVER MANUAL MEDIUM

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IMPLANT DRIVER SET ONE INTERNAL

The set includes the drivers suitable for the internal hexagon implant connection, available in mechanical and manual versions. For all implant driver variants, refer to the ONE INTERNAL catalog section.

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
GDS10

ONE CONICAL IMPLANT DRIVER SET

07GDD01	IMPLANT HANDPIECE DRIVER SHORT
07GDD02	IMPLANT HANDPIECE DRIVER LONG
02GDM10	IMPLANT MANUAL DRIVER

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IMPLANT DRIVER SET ONE CONICAL

The set includes the drivers suitable for the internal hexagon implant connection, available in mechanical and manual versions. For all implant driver variants, refer to the ONE CONICAL catalog section.

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ONE STEP DRILLS

The alveolar bone drills have a single color code ring to easily recognize the diameter. The alveolar bone drills have a smaller apex diameter to create the correct implant preparation with a soft and cancellous bone. Milled and laser marked work depth references. External irrigation drills. **Diamond Like Carbon Surface** finish.

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


OPTIONAL DRILLS




The optional drills allow to mark the bone for its preparation (marking burr) and easily drill the cortical bone (pointed starter drill). The Lindemann drill allows to correct the axis of the implant site after using the first surgical bone drill. The drill extension allows to increase the length of the handpiece tools.

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DRIVER PROTESICA**PROSTHETIC DRIVER****DESTORNILLADORES PROTESICO**

			
02EGM10 PROSTHETIC MANUAL DRIVER L = 10.0mm	07-EG05 PROSTHETIC HANDPIECE DRIVER SHORT L = 5.0mm	07-EG10 PROSTHETIC HANDPIECE DRIVER MEDIUM L = 10.0mm	07-EG15 PROSTHETIC HANDPIECE DRIVER LONG L = 15.0mm

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**PROSTHETIC DRIVER**

The **ONE implant system** has only one type of prosthetic driver (hexagon 1,2 mm) to screw and unscrew prosthetic components such as surgical screws, healing abutments, retaining screws for prosthetic components and transfer abutments. The prosthetic drivers feature the W&H hexagonal system.

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
MUA, OVERDENTURE DRIVER

Specific drivers for straight MUA abutments and overdenture abutments. The MUA drivers feature the W&H hexagonal system. The Uni-Base driver, used with the Uni-Base titanium bases and their specific retaining screws, allows to angle the insertion axis of the retaining screws up to 15° respect the implant axis.

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
**ACCESSORIES**

The **ONE system** features accessories to make easy specific surgical or prosthodontics operations. Parallel Pins are graduated and therefore can also be used to measure the osteotomy depth. The Depth Probe is graduated and has an apical cling to verify the position of the basal cortex in case of maxilla great sinus lift. The Torque Wrench has an adjustable screw in/out torque control (10-70 Ncm) and can be used both in the surgical and prosthodontics phase. The Ratchet Adapter allows you to manually use all handpiece tools and can be used within the torque wrench. If needed, Titanium Forceps allow you to handle the implant without any contamination. The Surgical Driver allows to position an implant or a retaining screw out of occlusion.

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BONE BOWL

The surgical kit can be completed with a bone bowl useful to correctly manage autologous bone or bone graft.

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
IMPLANT SIZING OVERLAY

The implant-sizing overlay allows the correct planning with TC, RX and OPT.

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ONE INTERNAL IS DESIGNED WITH THE MOST WIDESPREAD STANDARD O FINGER HEXAGON CONNECTION.

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A single prosthetic connection for every implant diameter.

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Machined platform and OsseoGRIP surface treatment up to 0,5 mm to implant platform to grant a safe juxta-crestal implant placement

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
Horizontal platform shifting to enhance the space available for peri-implant soft tissues

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


The implant is packaged without pre-assembled mounting device.
The standard cover screw is included in the package.
The self-centering design of cover screw allows easy insertion and tightening.
Use prosthetic driver to tighten the cover screw.
For more information on the implant, packaging and surgical procedures refer to the ONE Surgical Manual.

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IMPLANT DRIVERS

DESTORNILLADORES DE IMPLANTE

				
07NDD01 ONE INTERNAL HANDPIECE DRIVER SHORT LMAX = 5 mm	07NDD02 ONE INTERNAL HANDPIECE DRIVER LONG LMAX = 11 mm	07NDM10 ONE INTERNAL MANUAL DRIVER SHORT LMAX = 5 mm	07NDM11 ONE INTERNAL MANUAL DRIVER MEDIUM LMAX = 11 mm	07NDM12 ONE INTERNAL MANUAL DRIVER LONG LMAX = 16 mm

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IMPLANT DRIVERS

Implant drivers have a retention system to engage the implant, pick it up, carry it, and finalize positioning. Implant drivers feature the W&H hexagonal system. They have a reference that allows to control radial position of the implant's hexagonal index (indicating the position of the vertex of the hexagon).

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

COVER SCREW

TORNILLO DE SIERRE

Use Prosthetic Driver.
Tightening torque 10 Ncm.

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COVER SCREW

Surgical cover screw. Sterile packaging. Made in Ti gr.5.

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
HEALING ABUTMENTS

PILARES DE CICATRIZACION

Use Prosthetic Driver.
Tightening torque 20 Ncm.

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HEALING ABUTMENTS

They allow proper mucosal channel forming. All profiles are slightly enlarged in their emerging diameter respect the corresponding prosthetic abutments to allow an easy positioning, comfortable for the patient without use of anesthetic.

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
IMPRESSION POST DOUBLE TECHNIQUE

Impression post to be used with open and closed technique. This component has a series of accessories that allow both uses:

- Plastic cap with retentive notches to be placed over the abutment in such a way that the locking indexes are perfectly coupled to allow the laboratory a precise repositioning (closed technique). Made of POM.
- Cannula to extend the head of the retaining screw: it allows the user to facilitate the unscrewing of the retaining screw if the elastomer covers the head; subsequently remove the retentive body of the transfer abutment including an octagonal index for correct repositioning in the laboratory (open technique). Made of PMMA.

The impression post abutment and the implant retaining screw are made of Ti gr.5.

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IMPLANT REPLICA

Implant analog replicating the implant connection. Made in Ti gr.5.

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
IMPLANT SCAN ABUTMENT

PILAR DE ESCANEO PARA IMPLANTE

Use Prosthetic Driver.
Tightening torque 20 Ncm.
Retaining screw TYPE 1 included.

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IMPLANT SCAN ABUTMENT

Abutment for digital impression. Made in Ti gr.5 sandblasted.

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

UNI-BASE ABUTMENT

PILAR UNI-BASE

Use Prosthetic Driver.
Tightening torque 30 Ncm.
Retaining screw TYPE 1 included.

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UNI-BASE ABUTMENT

Uni-Base is the titanium base for CAD/CAM rehabilitation. It contributes to low cost workflows for the production of custom abutments.

- Accurate digital acquisition of the implant position through the scan abutment: extraoral on the model or intraoral.
- Designing abutment with CAD software.
- Realization of abutment mesostructure or crowns with CAM milling machine.
- Adhesive connection between Uni-Base and mesostructure or crown abutment.
- Available in different heights to allow maximum flexibility for customized restoration and with respect for biology.

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

TITANIUM ABUTMENT

PILAR EN TITANIO


Use Prosthetic Driver.
Tightening torque 30 Ncm.
Retaining screw TYPE 2 included.


01NMP10

TEMPORARY TITANIUM ABUTMENT
LOCKNG
0°
H = 1.7 mm


01NMP20

TEMPORARY TITANIUM ABUTMENT
NOT LOCKNG
0°
H = 1.7 mm

Use Prosthetic Driver.
Tightening torque 30 Ncm.
Retaining screw TYPE 2 included.
Angled respect hex vertex.

07MFN08

TITANIUM ABUTMENT
NOT LOCKNG
0°
SLEEVE H = 8,5 mm

07MFN12

TITANIUM ABUTMENT
NOT LOCKNG
0°
SLEEVE H = 12,5 mm

15MFN08

TITANIUM ABUTMENT
NOT LOCKNG
15°
SLEEVE H = 8,5 mm

15MFN11

TITANIUM ABUTMENT
NOT LOCKNG
15°
SLEEVE H = 11,5 mm

25MFN08

TITANIUM ABUTMENT
NOT LOCKNG
25°
SLEEVE H = 8,5 mm

25MFN11

TITANIUM ABUTMENT
NOT LOCKNG
25°
SLEEVE H = 11,5 mm

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TITANIUM ABUTMENT


Pre-finished abutments that allow fast prosthesis with high aesthetic value. The pre-finished angle of the abutment allows to easily correct the disparallelism of the implant. The sleeve of the abutments has an anti-rotational flat face and an anatomical profile to reduce the space occupied. Made of Ti gr.5.

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05MNA33


MULTI UNIT ABUTMENT
LOCKNG
H = 3,5 mm
30°

Use MUA Driver.
Tightening torque 30 Ncm.
Carrier included.

Use Prosthetic Driver.
Tightening torque 30 Ncm.
Retaining screw TYPE 3 included.
Angled respect hexagon vertex.
Carrier included.

Use Prosthetic Driver.
Tightening torque 30 Ncm.
Retaining screw TYPE 4 included.
Angled respect hexagon vertex.
Carrier included.

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MUA

Multiprotesic system with screwed connection. For screwed bridges, screwed retentive bars, Toronto-Bridge, totally edentulous arch rehabs. Ideal for rehabilitation with more than 3.0mm of gingiva. Compatible with the most popular standard on the market. Available straight and pre angled at 17° and 30°. Packed with carrier. Sterile packaging. Made in Ti gr.5.

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
COMPONENTES PROTESICOS MUA

01CMU10
MUA CASTABLE SLEEVE


Made in PMMA.
Use Prosthetic Driver.
Tightening torque 15 Ncm.
Retaining screw TYPE D included.

01TMU10
MUA IMPRESSION POST

01AMU10
MUA ANALOG


Made in Ti gr 5.
Use Prosthetic Driver.
Tightening torque 7 Ncm.

01CMU20
MUA TITANIUM SLEEVE

01GMU10
MUA HEALING CUP

Made in Ti gr 5.
Use Prosthetic Driver.
Tightening torque 15 Ncm.
Retaining screw TYPE D included.

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MUA PROSTHETIC COMPONENTS

The MUA system has a series of prosthetic components that allow to condition soft tissues, make an impression with open tray technique, create a master model with a dedicated analog and realize temporary and definitive prosthetic rehabilitation.

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

PRE-MILLED ABUTMENT

This prosthetic component is equipped with a premilled connection with all the characteristics and all the strict tolerances of implant system guarantee by ADVAN. The Pre-Milled abutment allows you to have customized transgingival design, regardless implant position; moreover this prosthetic system reduces production times and allows to orient each insertion axis of the secondary prosthesis.

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PROSTHETIC SCREWS

Prosthetic screw for components. Made in Ti gr.5.

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ONE CONICAL IS DESIGNED WITH THE GTB CONICAL COUPLING CONNECTION TO PERFORM EVEN WHERE TRADITIONAL CONNECTIONS STOP

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
A single prosthetic connection for every implant diameter

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Machined platform and OsseoGRIP surface treatment up to 0,2 mm to implant platform to grant a safe juxta-crestal implant placement

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

The implant is packaged without pre-assembled mounting device.
 The cover screw is included in the package.
 The self-centering design of cover screw allows easy insertion and tightening.
 Use Prosthetic Driver to tighten the cover screw.
For more information on the implant, packaging and surgical procedures refer to the ONE Surgical Manual.

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

DESTORNILLADORES DE IMPLANTE

02GDM10
IMPLANT MANUAL DRIVER
L_{MAX} = 9.0mm

07GDD01
IMPLANT HANDPIECE DRIVER
SHORT
L_{MAX} = 9.0mm

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IMPLANT DRIVERS

Implant drivers have a retention system to engage the implant, pick it up, carry it, and finalize positioning. Implant drivers feature the W&H hexagonal system. Implant drivers have laser marked references that indicate the transmucosal heights of pre-defined prosthetic abutments (positioned at 2.0-3.5-4.5-5.5 mm). They also have a reference that allows to control radial position of the implant's octagonal index (indicating the position of the vertex of the octagon).

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



HEALING ABUTMENTS

PILARES DE CICATRIZACION

Use Prosthetic Driver.
Tightening torque 7 Ncm.

		
050MG12 HEALING ABUTMENT H = 2.0 mm	050MG13 HEALING ABUTMENT H = 3.0 mm	050MG14 HEALING ABUTMENT H = 4.0 mm

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
IMPLANT SCAN ABUTMENT

PILAR DE ESCANEO PARA IMPLANTE

Use Prosthetic Driver.
Tightening torque 7 Ncm.
Retaining screw TYPE A included.

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IMPLANT SCAN ABUTMENT

Abutment for digital impression. Made in Ti
gr.5 sandblasted.

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

UNI-BASE ABUTMENT

PILAR UNI-BASE

Use Prosthetic Driver.
Tightening torque 25 Ncm.
Retaining screw TYPE A included.

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
TITANIUM ABUTMENT

PILAR EN TITANIO

Use Prosthetic Driver.
Tightening torque 25 Ncm.
Retaining screw TYPE A included.
Angled respect octagon vertex.

 07MFO08 TITANIUM ABUTMENT NOT LOCKNG 0° SLEEVE H = 8,5 mm	 15MFO08 TITANIUM ABUTMENT NOT LOCKNG 15° SLEEVE H = 8,5 mm	 25MFO08 TITANIUM ABUTMENT NOT LOCKNG 25° SLEEVE H = 8,5 mm
07MFO12 TITANIUM ABUTMENT NOT LOCKNG 0° SLEEVE H = 12,5 mm	15MFO11 TITANIUM ABUTMENT NOT LOCKNG 15° SLEEVE H = 11,5 mm	25MFO11 TITANIUM ABUTMENT NOT LOCKNG 25° SLEEVE H = 11,5 mm

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
H = 3,5 mm

30°

Use MUA Driver.
Tightening torque 25Ncm.
Mount pre-assembled included.

Use Prosthetic Driver.
Tightening torque 25 Ncm.
Retaining screw TYPE A included.
Angled respect octagon vertex.
Mount pre-assembled included.

Use Prosthetic Driver.
Tightening torque 25 Ncm.
Retaining screw TYPE C included.
Angled respect octagon vertex.
Mount pre-assembled included.

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PROSTHETIC SCREWS

Prosthetic screw for components.
Made in Ti gr.5.

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The aim of **Adván** is to deliver to the treatment Providers products that guarantee safety, efficacy and allow professional efficiency. During our activities for product development and product improvement **Adván** reserves the right to make changes to the products, instruments, surgical techniques and Instructions for use (IFU) if it constitutes a valuable improvement for the treatment Providers. **Adván** reserves the right to modify prices, delivery terms and sales terms without notice. Products might not be available in all Countries.

Lifetime Warranty: all implants and prosthetic components of **Adván** include a lifetime warranty that guarantees to replace any product that fails both for materials or manufacturing defects. This warranty is valid only if products are used in conjunction with the associated original components and instruments, used according to the Instructions for Use (IFU). The warranties given are exclusive for the benefit of eligible treatment Providers and are not for any other person or entity, including any patient.

Adván repairs or replaces the defective products according to the product defect and corrective action requested. The treatment Provider assumes all risks and liabilities deriving from the use of product, the lifetime warranty is limited to the replacement of defective product and **Adván** disclaims any liability to the treatment Provider or Patient and does not compensate costs for any additional components, instruments, treatment costs and any other cost related to or originated from the replacement of defective product.

Treatment planning and clinical application of the **Adván** products are under exclusive responsibility of treatment Provider; **Adván** recommends to attend intensive training programs, to complete a postgraduate dental implant courses and to strict adherence the Instructions for Use (IFU) and any other indication contained in the technical documentation available to Users. The product sale is authorized only to treatment Providers according with local regulatory rules and Laws. The products of this catalog could be patented or patent pending. Any replacement of product can be authorized by **Adván** only after submission of the "Guarantee Questionnaire" available upon request from Customer Service Office.

The "Guarantee Questionnaire" form must be filled in all the parts, duly signed and stamped, accompanied by the failed product and the other components placed by the treatment Provider within ONE (1) month after the product failure occurred. In case of an implant failure, x-ray of the failed implant IS MANDATORY. Prior to submit implant and the other components used in the treatment, ensure that sterilization of all the product has been made.

Warning – read carefully instructions for use (IFU) before product use

LIMITS AND LIMITATIONS

ADVÁN HEREBY DISCLAIMS ANY OTHER WARRANTIES, EXPRESS OR IMPLIED AND **ADVÁN** HEREBY EXCLUDES ANY LIABILITY FOR LOST EARNINGS AND DIRECT OR INDIRECT DAMAGES AS WELL AS COLLATERAL AND CONSEQUENTIAL DAMAGES, DIRECTLY OR INDIRECTLY RELATED TO **ADVÁN** PRODUCTS, SERVICES OR INFORMATION.

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SYMBOLS DESCRIPTION

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	Lot number
	Reference Number
	Sterilized using irradiation
	Do not re-use
	Non sterile
	Use by-date
	Advan products covered by the CE mark fulfill the requirements of the Directive 93/42/CEE concerning medical devices and falls within Classes Im, Is, IIa, IIb
	Advan products covered by the CE mark without the identification number fulfill the requirements of the Directive 93/42/CEE concerning medical devices and falls within Class I
	Consult instruction for use
	Manufacturer
	Distributor
	Not locking prosthetic component
	Octagon locking prosthetic component
	Hex locking prosthetic component
	Multi packaging (the number reported in the symbol refers to the number of units in the packaging)

Advan s.r.l.

Via Fratelli Solari, 5
33020 Amaro (UD) - Italia
Tel. +39 0433.096245
Fax +39 0433.096246

info@advanimplantology.com

www.advanimplantology.com